




Results

DO NOT COMPARE RESULTS WITH OTHER THERMOMETERS

Please interpret the results using this specific chart:

	34-37.1°C (93.2-98.8°F)	NORMAL
	37.2-38.1°C (98.9-100.5°F)	HIGH
	38.2-43°C (100.6-109.4°F)	FEVER

Accuracy

IT IS NORMAL THAT FOREHEAD THERMOMETERS GIVE DIFFERENT RESULTS THAN TRADITIONAL ONES




WARNING

Different parts of the body have different temperatures, so temperature readings taken from the forehead, mouth, armpit, ear, or anus will vary and should not be directly compared to each other, even if they are taken at the same time.

Résultats

**NE PAS COMPARER LES RÉSULTATS
AVEC D'AUTRES THERMOMÈTRES**

**Veillez interpréter les résultats à l'aide de ce
tableau:**

	34-37.1°C (93.2-98.8°F)	NORMALE
	37.2-38.1°C (98.9-100.5°F)	HAUTE
	38.2-43°C (100.6-109.4°F)	FIÈVRE

Précision

**IL EST NORMAL QUE LES THERMOMÈTRES
FRONTAUX DONNENT DES RÉSULTATS
DIFFÉRENTS DES THERMOMÈTRES
TRADITIONNELS**

WARNING

Les différentes parties du corps ont des températures différentes, c'est pourquoi les mesures de température prises sur le front, la bouche, l'aisselle, l'oreille ou l'anus varient et ne doivent pas être directement comparées entre elles, même si elles sont prises en même temps.